

## CLINICAL APPLICATION OF HA+HD



Blood Purif. Published online August 11, 2022. doi:10.1159/000525225

### Randomized Control Study on Hemoperfusion Combined with Hemodialysis versus Standard Hemodialysis: Effects on Middle-Molecular-Weight Toxins and Uremic Pruritus<sup>[4]</sup>

438 MHD patients from 37 HD centers

*"The long-term HP(HA)+HD can reduce  $\beta$ 2M and PTH levels and improve pruritus in MHD patients independently on the use of high-or low-flux dialyzers, showing that the results are linked to the effect of adsorption."*

All results are based on specific patients, actual outcomes may vary individually.



Blood Purif (2022)51(3):199-212. https://doi.org/10.1159/000514940

### Survival Outcomes of Hemoperfusion and Hemodialysis versus Hemodialysis in Patients with End-Stage Renal Disease: A Systematic Review and Meta-Analysis<sup>[5]</sup>

*"HP(HA)+HD can achieve the complementary elimination of metabolites, preventing and treating complications caused by long-term dialysis."*

All results are based on specific patients, actual outcomes may vary individually.



2021 Jul;9(14):1133. doi:10.21037/atm-21-1100.

### Cost-effectiveness analysis of hemodialysis plus hemoperfusion versus hemodialysis alone in adult patients with end-stage renal disease in China<sup>[6]</sup>

*"Our analyses indicate a potential for HD+HP(HA) to be cost-effective for adult patients with ESRD in China. Compared to HD alone, HD+HP reduces incidence of severe CVD events."*

All results are based on specific patients, actual outcomes may vary individually.

For Professional Academic Exchange Only. Not for commercial use. Reproduction or distribution is strictly prohibited.



## PRODUCT SPECIFICATIONS

Adsorbent Volume (mL)	130
Biocompatibility	Tested as required in ISO10993
Adsorbent Material	Double Cross-linked Styrene-divinylbenzene Copolymers
Housing Material	Polycarbonate
Sterilization Method	Irradiation Sterilization

Precaution: Physicians should pay attention to potential side effects such as low blood pressure, chill and coagulation when using the product. \*Contraindications, Warnings and Precautions: refer to Instructions for Use.


## Jafron: Global Manufacturer and Supplier of Adsorption Columns



CE-MDR



ISO 9001



EN ISO 13485

### References

- [1] Cozzolino M, Magagnoli L, Ciceri P. Toxins (Basel). 2025 Jun 10;17(6):295
  - [2] Ronco Claudio. Integrative Medicine in Nephrology and Andrology 11(1):e00006, March 2024.
  - [3] Expertconsensus on the clinical application of hemoperfusion in maintenance hemodialysis patients(2025 Edition)[J]. SMJ, 2025, 48(9):529-542.
  - [4] Zhao D, et al. Blood Purif. 2022 Aug 11;1-11. doi: 10.1159/000525225. Epub ahead of print. PMID: 35952629.
  - [5] Cheng W, Luo Y, Wang H, Qin X, Liu X, Fu Y, Ronco C. Blood Purif. 2022;51(3):213-225. doi: 10.1159/000514187. Epub 2021 May 10. PMID: 33971651; PMCID: PMC8985078.
  - [6] Wang H, Jin H, Cheng W, Qin X, Luo Y, Liu X, Fu Y, Jiang G, Lu W, Jin C, Pennington M. Ann Transl Med. 2021 Jul;9(14):1133. doi: 10.21037/atm-21-1100. PMID: 34430574; PMCID: PMC8390641.
- For Professional Academic Exchange Only. Not for commercial use. Reproduction or distribution is strictly prohibited.

## JAFRON BIOMEDICAL CO., LTD.

Address: No. 98, Technology Sixth Road, High-tech Zone, Zhuhai City, 519085, Guangdong, China.  
Tel: +86 (756) 3689708  
E-mail: customerservice@jafron.com  
Website: www.jafron.com



## This product is imported and distributed by: Husk Medical B.V.

Address: Kelwinring 44, 2952 BG Alblasserdam, The Netherlands.  
Tel: +31 78 30 32 867  
Email: info@huskmedical.com  
Website: www.huskmedical.com



Stock Abbreviation: JFSW  
Stock Code: 300529

# HEMOADSORPTION TECHNOLOGY



## HA130

Order code Husk Medical B.V.: HM-675100

20+ YEARS IN BLOOD PURIFICATION CLINICAL PRACTICES  
USED IN MORE THAN 90 COUNTRIES AND DISTRICTS

For Professional Academic Exchange Only. Not for commercial use. Reproduction or distribution is strictly prohibited.

CY-HA130G02-02-202601-EN

# HEMOADSORPTION TECHNOLOGY

HA130 adsorbs uremic toxins accumulated in ESRD<sup>[1]</sup>

## Jafron Resin Technology


- Pore size adjustment**  
Extensive sieving of endogenous and exogenous molecules
- Hyper double cross-linked**  
High mechanical strength with surface area over 1000m<sup>2</sup>/g
- Nanobiomaterial coating**  
Good biocompatibility, preventing particle detachment into human body

Large-middle molecule toxins<sup>[1]</sup>

Medium-middle molecule toxins<sup>[1]</sup>  
TNF- $\alpha$ , IL-6, IL-1 $\beta$

Small-middle molecule toxins<sup>[1]</sup>  
iPTH, IL-8,  $\beta$ 2-MG

Protein-bound uremic toxins (PBUTs)<sup>[1]</sup>  
IS, PCS, Hcy



## Adsorption Mechanism

### 3D mesh molecular sieve

The resin can sieve uremic toxins, optimize removal of uremic middle molecular toxins and protein bound toxins in the blood<sup>[1]</sup>

### Lipophilic and hydrophobic effects

The resin surface contains abundant non-polar structures (such as benzene rings), which interact hydrophobically, with the nonpolar regions of toxins.


### Van der Waals force

Toxins are adsorbed onto an over-1000 m<sup>2</sup>/g surface in the porous micro-environment of the Jafron HA130 resin.

# HA+HD: COMBINED ARTIFICIAL KIDNEY

HA+HD Therapy Operation Mode

HA130 cartridge is compatible with many extracorporeal blood circuits including Hemodialysis (HD), Hemofiltration (HF) & Hemodiafiltration(HDF). Indicating accumulated toxins in end stage renal disease hemodialysis related complications.




## Comparison between different dialysis modalities<sup>[1]</sup>

Molecular Weight	Exogenous Molecules			Endogenous Water-Soluble Molecules		
	Gut-Derived Protein-Bound Molecules<80%	Gut-Derived, Protein-Bound Molecules>80%	Small Molecules (<0.5 kDa)	Small-Middle Molecules (0.5-15 kDa)	Medium-Middle Molecules (15-25 kDa)	Large-Middle Molecules (25-58 kDa)
Uremic toxins	ADMA, SDMA, uric acid, carbamylated compounds, urea, TMAO	Hcy, IS, pCS, CML, kynurenes	Urea	$\beta$ 2-microglobulin, IL-8	TNF, IL-18, IL-10, IL-6, kappa-FLC, myoglobin, sTNFR2, FGF-2, prolactin, complement factor D	AGEs, FGF-23, $\lambda$ -FLC
Low-flux HD	YES	NO	YES	NO	NO	NO
High-flux HD	YES	NO	YES	YES	NO	NO
Online HDF/HDX	YES	YES	YES	YES	YES	NO
Hemoadsorption (HA)	YES	YES	—	YES	YES	YES

HA+HD: Inflammation Mitigation Hypothesis<sup>[2]</sup>

The enhanced clearance of toxins and mediators leads to an increased elimination of  $\beta$ 2M.

3-Sessions-per-week regimen achieves a progressive reduction in pre-dialysis  $\beta$ 2M values.



# INDICATION OF HA+HD

Expert Consensus on The Clinical Application of Hemoperfusion in MHD Patients (2025 Edition)<sup>[3]</sup>

## 1 Standard Treatment Regimen

1 session per week

Complications	Appropriate Patients	Recommended frequency
Uremic Pruritis	Duo's pruritus score of the MHD patient >12 or the VAS score $\geq$ 7	Once weekly
Uremia-Related Sleep Disorders	Pittsburgh sleep quality index(PSQI)of the MHD patient $\geq$ 10	Once weekly / every 2 weeks
Restless Leg Syndrome(RLS)	RLS severity scale scores $\geq$ 11	Once weekly
Microinflammatory Condition	High-sensitivity C-reactive protein (hs-CRP) >3mg /L,CRP >8 mg /L, IL-6 $\geq$ 16pg/mL, or TNF-a $\geq$ 41 pg/mL	Once weekly
CKD-Mineral and Bone Disorder(CKD-MBD)	For MHD patients with moderate to severe SHPT uncontrolled by drugs, showing iPTH > 600 pg/mL or significant vascular calcification (Agatston >100 or Kauppila >4.5)	Once weekly
High serum $\beta$ 2-MG	$\beta$ 2-microglobulin( $\beta$ 2-MG)>30 mg/L or with dialysis related amyloidosis such as carpal tunnel syndrome	1-3x weekly
Protein-Energy Wasting (PEW)	The modified quantitative subjective global assessment(MQSGA)score >20 or the malnutrition inflammation score(MIS)>18	Once weekly
Refractory Hypertension	Systolic blood pressure(SBP)>160 mm Hg(1 mmHg=0.133kPa)	Once weekly / every 2 weeks
Uremic Peripheral Neuropathy	Limb numbness or paresthesia,weakened or absent muscle tension or tendon reflex or whose peripheral nerve electrophysiological examination shows more than 2 nerves	Once weekly
Protein-Bound Uremic Toxins(PBUTs) Accumulation	For MHD patients exhibiting accumulation of protein-bound uremic toxins (PBUTs), especially when indoxyl sulfate (IS) > 44 mg/L or p-cresyl sulfate (PCS) > 19 mg /L	Once weekly

## 2 Intensive Treatment Regimen

2-3 sessions per week

For Maintenance Hemodialysis (MHD) patients presenting with severe, refractory complications that are unresponsive to conventional HD/HDF/HDX, an individualized, short-term intensive regimen is recommended.