

flexGRIP®

catheter cushion and **stabilization device**

ALPHA

the answer to known issues
with sutures:

- ✓ Central Line Associated Blood Stream Infections (CLABSI) are reduced by 80% when using a fixation device ¹
- ✓ designed to improve patient comfort & safety
- ✓ no more needle stick injuries for care givers

safe
comfortable
easy to use
versatile

Find out more on www.flex-grip.com




bedal International

flexGRIP®

catheter cushion and **stabilization device**

The FlexGRIP Alpha is a single use patch for the stabilization of central venous catheters (CVC), peripherally inserted catheters (PICC), arterial lines and midlines containing wings with suture holes.



ALPHA

features

- hydrocolloid tape for sensitive skin (CVC / Arterial / Hemodialysis version)
- hypo-allergenic acrylic tape (PICC / Midline version)
- "suture-like" flexible posts can be closed/opened
- tape design enables use of Biopatch

benefits

- 2 posts replace the need for suturing
- compatible with most catheters with suture wings: CVC / PICC / Midline / Arterial / Hemodialysis)
- intuitive, fast and easy to use
- designed to prevent MARS!

instruction video



ordering information

ref. no.	order code Husk Medical	description & size range	qty
ALPHA00	HM-522078	Alpha PICC / Midline	50 pcs/box
ALPHASS	HM-522079	Alpha CVC / Arterial	50 pcs/box

M HUSK MEDICAL

This product is imported and distributed by
Husk Medical B.V.

Kelvinring 44
2952 BG Alblasterdam
Tel: 078-3032867

Email: info@huskmedical.com
www.huskmedical.com



¹ Crnich CJ, Maki DG. The promise of novel technology for the prevention of intravascular device-related bloodstream infection. II. Long-term devices. Clin Infect Dis. 2002; 34(10):1362-1368.