

flexGRIP®

catheter cushion and **stabilization device**

UNIVERSAL

The answer to known issues
with sutures:

Designed to:

- ✓ minimize complications associated with suturing
- ✓ improve patient comfort and safety
- ✓ no more needle stick injuries for care givers

safe
comfortable
easy to use
versatile

Find out more on www.flex-grip.com



HIGH FIXATION STRENGTH
WITH
AIR PILLOW
TECHNOLOGY
NO IMPACT TO FLOW RATE


bedal International

flexGRIP®

catheter cushion and **stabilization device**

suitable for stabilizing & fixating a large range of catheters, such as:

nephrostomy (post-op),
drainage catheters, gripper
needles (Port-a-cath),
thrombolysis catheters,
jejunostomy, tubings, cables
and drivelines (eg. LVAD).



UNIVERSAL

features

- large diameter range: 5-18 Fr
- unmatched **strength** of securement
- **easy to use**; pre-mounted fixation leash
- additional lock for leash
- **skin-safe** hypoallergenic adhesive
- **soft** materials only

benefits

- universal system; wide variety of applications
- designed to:
 - prevent catheter migration & related complications
 - prevent MARSII – medical adhesive related skin injury
 - optimize patient comfort
- extra safety lip prevents accidental loosening
- extra patient comfort

instruction video



ordering information

ref. no.	order code Husk Medical	description & size range	qty
UN511	HM-522075-S	FlexGRIP Universal 5-11 Fr	50 pcs/box
UN1118	HM-522075-M	FlexGRIP Universal 11-18 Fr	50 pcs/box

M HUSK MEDICAL

This product is imported and distributed by
Husk Medical B.V.

Kelvinring 44
2952 BG Alblasterdam
Tel: 078-3032867

Email: info@huskmedical.com
www.huskmedical.com



bedal International