


Protective covers for positioning pillow and birthing ball

Single-use covers for hygiene during childbirth



REF 808100



REF 808101

Product description

These single-use covers provide a hygienic barrier between the user and a positioning pillow or birthing ball.

The positioning pillow cover (REF 808100) features a zipper for easy fitting and removal. The birthing ball cover (REF 808101) features an elasticated opening for easy fitting and removal and to help keep the cover securely in place during use.

Both covers are made from blue SMS nonwoven material and are supplied non-sterile. Replace the cover after each user.

Product specifications

Material	SMS (Spunbond Meltblown Spunbond)
Material weight	35 g/m ² (+/- 3 g/m ²)
Sterile	No
Single-use	Yes
Colour	Blue
Intended use	REF 808100: Single-use protective cover for a positioning pillow, intended to provide a hygienic barrier between the user and the pillow. REF 808101: Single-use protective cover for a birthing ball, intended to provide a hygienic barrier between the user and the ball.
Closure	REF 808100: Zipper REF 808101: Elasticated opening
Applicable norms	• EN ISO 15223-1:2016 • EN ISO 14971:2019



Quality and certification

Classification	Medical Device class I
Certifications	• Manufactured as per ISO 13485:2016 • EU Declaration of Conformity (MDR)
Notified Body (CE)	Not applicable for non-sterile class I medical devices without a measuring function

Ordering information

REF code	Product description	Packing unit
808100	Protextra™ cushion cover type 'Boomerang', SMS material, with zipper, blue color	BOX/75
808101	Protextra™ protective cover for birthing ball, SMS material, with elastic, blue color	BOX/75

Document code: PDS/EN/PPC/1.0

Version: 1.0

Last update: 06-11-2025

PROTEXTRA™ is a registered brand of Husk Healthcare Group

Husk Medical B.V.

Kelvinring 44

2952 BG Alblasserdam,

The Netherlands

Mail: info@huskmedical.com