

SteriFlow™

SPINAL NEEDLES (QUINCKE & PENCIL POINT)

Precision and convenience during spinal anaesthesia and lumbar puncture




Product description

Precision is essential during subarachnoid anesthesia and diagnostic lumbar punctures. The needle must provide reliable access to the subarachnoid space with minimal tissue trauma and attention to patient comfort. A stable, well-controlled puncture supports the healthcare professional in performing the procedure efficiently and carefully.

SteriFlow™ spinal needles have been developed for the precise execution of spinal anesthesia and diagnostic punctures. They are available in both Quincke and Pencil Point versions, allowing the choice to be tailored to clinical preference, puncture technique, and patient situation.

The Quincke needle has a sharp, cutting tip and is suitable for a direct and controlled puncture. The Pencil Point needle is designed with an atraumatic tip, aimed at spreading rather than cutting the dural fibres. This can help limit tissue trauma and patient discomfort after the puncture. Scientific literature shows that atraumatic Pencil Point needles are associated with a lower incidence of post-dural puncture headache compared to conventional cutting needles. Thus, the Pencil Point version supports a patient-centered approach in procedures where precision, control, and comfort are important.¹



Order information

Quincke tip without introducer



Product code	Diameter	Length (inch/mm)	Color coding
400510	16G	3½" / 90mm	Light gray
400511	18G	3½" / 90mm	Pink
400512	19G	3½" / 90mm	Light brown
400513	20G	3½" / 90mm	Yellow
400514	21G	3½" / 90mm	Green
400515	22G	3½" / 90mm	Dark gray
400516	23G	3½" / 90mm	Dark blue
400517	24G	3½" / 90mm	Purple
400518	25G	3½" / 90mm	Orange
400519	26G	3½" / 90mm	Brown
400520	27G	3½" / 90mm	Light blue

Quincke tip with introducer



Product code	Needle diameter	Needle length (inch/mm)	Introducer diameter	Introducer length (inch/mm)	Color coding
400521	22G	3½" / 90mm	19G	1½" / 38mm	Dark gray
400522	23G	3½" / 90mm	19G	1½" / 38mm	Dark blue
400523	24G	3½" / 90mm	20G	1½" / 38mm	Purple
400524	25G	3½" / 90mm	20G	1½" / 38mm	Orange
400525	26G	3½" / 90mm	21G	1½" / 38mm	Brown
400526	27G	3½" / 90mm	22G	1½" / 38mm	Light blue

Pencil Point tip without introducer



Product code	Diameter	Length (inch/mm)	Color coding
400530	18G	3½" / 90mm	Pink
400531	19G	3½" / 90mm	Light brown
400532	20G	3½" / 90mm	Yellow
400533	21G	3½" / 90mm	Green
400534	22G	3½" / 90mm	Dark gray
400535	23G	3½" / 90mm	Dark blue
400536	24G	3½" / 90mm	Purple
400537	25G	3½" / 90mm	Orange
400538	26G	3½" / 90mm	Brown
400539	27G	3½" / 90mm	Light blue



Pencil Point tip with introducer



Product code	Needle diameter	Needle length (inch/mm)	Introducer diameter	Introducer length (inch/mm)	Color coding
400540	22G	3½" / 90mm	19G	1½" / 38mm	Dark gray
400541	23G	3½" / 90mm	19G	1½" / 38mm	Dark blue
400542	24G	3½" / 90mm	20G	1½" / 38mm	Purple
400543	25G	3½" / 90mm	20G	1½" / 38mm	Orange
400544	26G	3½" / 90mm	21G	1½" / 38mm	Brown
400545	27G	3½" / 90mm	22G	1½" / 38mm	Light blue

Packaging information

Sterile individually packed in a blister pack, then per BOX/50.

Scan this QR code to visit the product page for more information:



SteriFlow™

STERIFLOW™ is a registered trademark of **Husk Healthcare Group B.V.**

Interested?

Contact your local partner for more information, by scanning the QR code and selecting your country:

