SteriFlow™

HYPODERMIC INJECTION NEEDLES


Very wide range of diameters and lengths



SteriFlow™

Product description

SteriFlow™ injection needles are made of high-quality steel and feature a semi-transparent hub to achieve optimal flow observation. The specially designed needle tip ensures optimal skin penetration with minimal impact and the lowest possible pain level. All needles are color-coded in accordance with ISO 6009:2016 for easy and standardized identification.



Simple color-coded identification

Optimal flow observation through semi-transparent hub

Low penetration resistance due to consistent optimal needle grinding

Specifications

Classification	Medical Device Class IIa	**Packaging information**	Sterile individually packed in a blister pack, then in BOX/100.
Certification	· Produced according to ISO 13485:2016 · CE certified (MDR EU 2017/745) · EU Declaration of Conformity (MDR EU 2017/745)	**Other specifications**	· Compatible with all standard Luer tip syringes · Latex-free, sterile, and for single use · Individually sterile packed in blister packaging

Ordering information

Scan the QR code to go to the product page with an up-to-date overview of the options.



SteriFlow™

STERIFLOWTM is a registered trademark of **Husk Healthcare Group B.V.**



Interested?
Contact your local partner for more information, by scanning the QR code and selecting your country: